# **Science Together**




# Mounting Bracket AZURA® Bio L

The mounting bracket can be mounted onto the right or left side of AZURA® devices with a L footprint. It can hold one manual injection valve, one pH flow cell and a column at the same time.


### Attaching a mounting bracket to the side panel

**Tool** • Allen screwdriver, size 2.5 mm



### Mounting a column holder

**Tool** • Allen screwdriver, size 2 mm



# Mounting a pH flow cell

**Tool** • Torx screwdriver, T10

Process	Figure
<ol> <li>Lay the pH flow cell against the upper inclined positioned holes of the mounting bracket.</li> <li>Take the two screws ① out of the pH flow cell accessories kit.</li> <li>Fasten the pH flow cell from behind with the two screws.</li> </ol>	
	Fig. 3: Mounted pH flow cell


## Mounting a valve

**Tool** • Allen wrench, size 2.5 mm

■ Torx screwdriver, T20

Process	Figure
<ol> <li>Remove the lid from the lever.</li> <li>Unscrew the screw of the lever ① with the screwdriver size 2.5 mm.</li> <li>Pull off the lever from the adapter plate ②.</li> </ol>	
<ul><li>4. Loosen the 2 screws of the adapter plate 2 with the torx screwdriver T20.</li><li>5. Detach the adapter plate 2 from the feather key 3.</li></ul>	
6. Detach feather key ③ from valve head ④.	3
	Fig. 4: Exploded view

# 7. Place the valve head on the front of the mounting bracket so that the mark ⑤ points upwards and the groove ⑥ (fig. 5) points to the left.



Process	Figure
<ul> <li>8. Mount the feather key ⑦ pointing to the AZURA® L system onto the groove ⑥.</li> <li>9. Mount the adapter plate ⑧ with the cable pointing down onto the feather key ⑦ and tighten it with the screwdriver T20.</li> <li>10. Tighten the lever with the screwdriver size 2.5 mm and mount the lid onto the lever.</li> </ul>	
	Fig. 6: Mounting a valve



**Note:** In load position, the 6-port valve is connected: port 1 with port 6, port 2 with port 3 and port 4 with port 5.