Smartline

► Manager 5050 Manual/Handbuch

V7599

Table of Contents

Note: For your own safety, carefully read the manual and always observe the warnings and safety information on the device and in the manual!

Using this Manual	7
Conventions in this manual	7
SOP's in this manual	7
Smartline Manager 5050	8
HPLC	8
Manager - Device types	8
Laboratory use	8
Where is it prohibited to use the device or system?	8
Safety	9
Laboratory regulations	9
Solvents	9
PEEK connections	11
Protective measures	11
Target group	12
Operating the device or device system	12
To what should the user pay particular attention?	12
What expertise should users have to safely operate an HPLC device or device system?	12
Symbols and labels	13
Installation	14
Packaging and transport	14
Fastening material and shipping boxes	14
Scope of supply	14
Device types	14
Accessories	14
Checking the scope of delivery	14
Space requirements	15
Installation site	15
Low Pressure Gradient (LPG) Block	16
General Description	16
Installation	16
Installing the manifold on the pump head	17
Electrical Connection with the Smartline Pump 1050	18
Eluent Connection with the Smartline Pump 1050	18
Vacuum Degasser Unit	20
Introduction	20
Working Principle	20
Principles of degassing using Teflon AF [®] membranes	20
General Description	21
Connection	21
Front Panel Connections	21
Front Panel LED's	22
Rear Panel Connections	23
System Requirements	23
Solvents/Mobile Phase	23
Chemical compatibility	23
Space Requirements	24
Electrical Power Requirements	24

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Using this Manual

Conventions in this manual


- 1 **Special Warnings** are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates **helpful advices**.

Important Hints are marked by the marginal hand symbol.

Arrows in an outlined form like this example used in block diagrams, indicate an automatic program run and change to the next line without the necessity of manual interventions.

Arrows like this, used in block diagrams, indicate that the user is asked to press the corresponding arrow keys

The triangles symbolize the use of corresponding arrow keys.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your Smartline Manager 5050. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and be copied, applied signed, and filed in order to document the performance of the instrument.

Please operate the instrument and all accessories according to instructions and SOP's in this manual. This ensures proper results and longevity of your equipment.

Table 1 *SOP's in this manual*

SOP 1	Selection of the Gradient Mode	18
SOP 2	Selection of operation mode:	19

Smartline Manager 5050

HPLC

High-performance liquid chromatography (HPLC) is a method for separating substance mixtures, determining substances and measuring their concentration.

The device or system is suitable for high-performance liquid chromatography. It is suitable for laboratory use, for analyzing substance mixtures that can be dissolved in a solvent or solvent mixture.

Manager - Device types

- Smartline Manager with degasser, low pressure gradient (LPG) and interface module, available also preparative or bio compatible.
- Smartline Manager with degasser, low pressure gradient module (LPG), available also bio-compatible.
- Smartline Manager with degasser module and interface module

Table 2 Smartline Manager 5050

Smartline Manager 5050					
Interface Module	Degasser Module Analytical	Degasser Module Preparative	LPG Module	Biocompatible	Order number
x	x		x		A5331
x	x		x	x	A5331-1
x	x				A5332
	x		x		A5333
	x		x	x	A5333-1
x		x	x		A5338
x		x	x	x	A5338-1

Laboratory use

- Biochemistry analyses
- Chemical analyses
- Food analyses
- Pharmaceutical analyses
- Environmental analyses

Where is it prohibited to use the device or system?

Risk of explosion! Never use the device in potentially explosive atmospheres without appropriate protective equipment and approval by a notified body! Inform the technical support department of the manufacturer!

Safety

Laboratory regulations

Observe national and international regulations pertaining to laboratory work!

- Good Laboratory Practice (GLP) of the American Food & Drug Administration
- For development of methods and validation of devices: Protocol for the Adoption of Analytical Methods in the Clinical Chemistry Laboratory, American Journal of Medical Technology, 44, 1, pages 30–37 (1978)
- Accident prevention regulations published by the accident insurance companies for laboratory work

Solvents

Solvents suitable for use in HPLC:

- Acetate buffer solutions
- Acetone at 4–25 °C (39.2–77.0 °F)¹
- Acetonitrile²
- Benzene
- Carbon dioxide (liquid 99.999 % CO₂)
- Chloroform
- Dilute acetic acid (e.g. 0.1-1 %) at 25 °C/77.0 °F
- Dilute ammonia solution
- Dilute sodium hydroxide (1 M)
- Ethyl acetate
- Ethanol
- Formate buffer solution
- Isopropanol
- Methanol
- Phosphate buffer solutions (0.5 M)
- Toluol
- Water

Even small quantities of other substances, such as additives, modifiers, or salts can influence the durability of the materials. The list of selected solvents was compiled based on research in the pertinent literature and is only a recommendation by the manufacturer. In the event of doubt, contact the technical support of the manufacturer.

To avoid damage from leaks, always place solvent bottles in a solvent tray on the device.

¹ valid for the specified temperature range

² not recommended in combination with PEEK small parts or PEEK capillary

Organic solvents are toxic above a certain concentration. Ensure that work areas are always well-ventilated! Wear protective gloves and safety glasses when working on the device!

Connect silicon tube (inner diameter: 3.0 mm) with the olive-type tube fitting of the exhaust and lead the gases or liquids into a suitable collecting container or to a fume hood.

Organic solvents are highly flammable. Since capillaries can detach from their screw fittings and allow solvent to escape, it is prohibited to have any open flames near the analytical system!

Only use solvents that have a self-ignition point higher than 150 °C under normal ambient conditions!

The following solvents can attack the components of the manager and are therefore not suitable:

- Concentrated mineral and organic acids
- Concentrated bases
- Eluents containing particles
- Halogenated hydrocarbons, e.g. Freon[®]
- Hexafluoroisopropanol
- Hexanes (60 % n-Hexane)
- Hydro fluoro eluents
- Perfluorinated eluents, e.g. Fluorinert[®] FC-75, FC-40
- Perfluorinated polyether, e.g. Fomblin[®]

The following solvents are suitable to only a limited extent for use in the manager:

- Diethylamine (0.1 %) (DEA)
- Dilute phosphoric acid
- Dimethyl sulfoxide (DMSO)
- Methylene chloride³
- Slightly volatile eluents
- Tetrahydrofuran (THF)
- Triethylamine (0.1 %) (TEA)
- Trifluoroacetic acid (0.1 %) (TFA)

The following solvents are not suitable for use in the **degasser module**:

- Azides
- Benzene
- Carbon dioxide (liquid 99.999 % CO₂)
- Concentrated mineral and organic acids
- Concentrated bases
- Dilute sodium hydroxide (1 M)
- Eluents containing particles
- Halogenated hydrocarbons, e.g. Freon[®]
- Hexafluoroisopropanol

³ not recommended in combination with PEEK small parts and PEEK capillary

- Hexanes (60 % n-Hexane)
- Hydro fluoro eluents
- Perfluorinated eluents, e.g. Fluorinert[®] FC-75, FC-40
- Perfluorinated polyether, e.g. Fomblin[®]

Regularly check for leaks and clogged capillaries – test back pressure without column!

PEEK connections

- Disposable PEEK fittings
- One-piece disposable polyetheretherketone fittings (PEEK) for easier installation of flexible, thin capillaries (Outer diameter 0.5 mm)
- Tightening torque of the PEEK screw:
Tightened by hand (approx. 0.5 Nm)

Protective measures

1. You are only permitted to perform the maintenance tasks described in this manual.
2. All other maintenance tasks are to be performed exclusively by the manufacturer or a company authorized by the manufacturer.

Without exception, the following applies to all maintenance tasks that can be performed by the user:

1. Switch off the device and pull the power plug!
2. Never open the device! High voltage poses a life-threatening risk!

Target group

Operating the device or device system

The device can be operated by means of the chromatography software at the workstation or the membrane keyboard at the device.

To what should the user pay particular attention?

To make your HPLC separations as efficient as possible, pay close attention to the following:

1. Once they have been used, never re-use capillaries in other areas of the HPLC system.
2. Only use a given fitting for one specific port and never re-use it for other ports. Always install new fittings on each separate port.
3. Use special columns and follow the manufacturer's instructions on caring for the columns!
4. Regularly check for clogged capillaries – test back pressure without column!
5. Use ultra-pure, filtered solvents – Gradient grade – for the HPLC system.
6. Filtration of substances to be analyzed.
7. Use of inline filters

Only allow the technical support department of the manufacturer or a company authorized by the manufacturer to open the devices for maintenance and repair work.

What expertise should users have to safely operate an HPLC device or device system?

- Completed degree as chemical laboratory technician or comparable vocational training
- Fundamental knowledge of liquid chromatography
- Participation in an installation of the system performed by the manufacturer or a company authorized by the manufacturer, or suitable training on the system and chromatography software
- Basic knowledge of Microsoft Windows®
- Knowledge regarding substances that are suitable only to a limited extent for use in liquid chromatography

Symbols and labels

Explanations of symbols and labels on the device or system

Hazard symbol indicating microelectronic devices that can be damaged by electrostatic discharge when touched.

Marking for devices that comply with the Canadian requirements for laboratory equipment: CAN/CSA-C22.2 No. 61010-1, second edition, including Amendment 1, or a later version

CE-marking for devices that comply with the pertinent EU directives (Conformité Européenne) and come with a declaration of conformity from the manufacturer.

For your own safety, read the manual and always observe the warnings and safety information on the device and in the manual!

Installation

Packaging and transport

At the factory, the device was carefully packed for safe transport.

Check the device for signs of damage that occurred during transport. If the shipment is incomplete or damaged, inform the manufacturer within three working days. Also inform the freight carrier about transport damage.

Fastening material and shipping boxes

The device is held in place and protected by foam inserts at the top and bottom. Please keep the transport box and the foam inserts.

- Remove the foam insert on the top of the device.
- Grip the device at its sides, near the front panel, and lift it out of the packaging.

Scope of supply

Device types

- Smartline Manager with degasser module, low pressure gradient module (LPG) and interface module
- Smartline Manager with degasser module and , low pressure gradient module (LPG)
- Smartline Manager with degasser module and interface module

Accessories

- Manual
- Kit for the manager

Only use original parts and accessories made by the manufacturer or a company authorized by the manufacturer.

Checking the scope of delivery

1. Check whether the device and accessories are complete.
2. If anything is missing, inform the technical support department of the manufacturer.

Hotline of the technical support department of the manufacturer:

Languages: German und English
telefonisch erreichbar: 8 – 17Uhr (CET)
Phone: +49-(0)30-809727-0
Telefax: +49(0)30-8015010
E-Mail: info@knauer.net (manufacturer)

Space requirements

- Side clearance to other devices:
- If there is a device on one side, min. clearance of 5 cm.
- If there are devices on both sides, min. clearance of 10 cm.

Make sure that the power plug on the rear of the device is always accessible, so that the device can be disconnected from the power supply

Installation site

- Air humidity: Below 90% (non-condensing)
- Temperature range: 4–40 °C; 39.2–104 °F
- Sunlight: When setting up the device at the installation location, make sure that it is protected against direct sunlight.

Low Pressure Gradient (LPG) Block

General Description

The LPG module in the Smartline Manager 5050 complements the Smartline Pump 1050 that creates a completely quaternary low pressure gradient system. The flow range reaches up to 50 ml/min depending on the pump head. The valve block is controlled over the Smartline Pump 1050 and delivers gradient mixtures with the highest precision over the entire flow range.


The LPG Block is integrated in the Smartline Manager 5050 which can be additionally equipped with a degasser and/or an interface module.

The LPG module in the Smartline Manager 5050 can only be operated in combination with a Smartline Pump 1050 or 1000. The Smartline Pump 1050 controls all electrical functions on the Smartline Manager 5050 and provides the required electrical power supply. This is realized with a connection cable (G0649). The connection cable is included in the accessories for the Smartline Manager 5050.

Control of the LPG block can be obtained in stand-alone mode over the Smartline Pump 1050 or with ChromGate[®] or ClarityChrome[®] software.


Installation

The Smartline Manager 5050 and the Smartline Pump 1050 must be arranged in a tower on top the other. See following figures:

- | | |
|-----------|--|
| 1.1 | Connection degasser to LPG unit |
| 1.2 | Static mixing chamber |
| 1.3 | Purging outlet |
| 1.4 – 1.7 | Connections LPG unit to pump head mixing block A-D |
| 1.8 | Manifold |

Figure 1 Arrangement of the Smartline Manager 5050 and Smartline Pump 1050 – front view


- 2.1 | Manager connector
- 2.2 | Connection cable
- 2.3 | Pump connector

Figure 2 Arrangement of the Smartline Manager 5050 and the Smartline Pump 1050 – back panel view

Installing the manifold on the pump head

- Connect up to four outlets of the manager LPG module to the intake manifold on the pump head.
- Close unused inlets by means of filler caps to prevent air from being sucked into the system.

- 3.1 | Knurled-head screw
- 3.2 | Filler cap

Figure 3 Intake manifold

1. Screw the intake manifold into the lower inlet fitting of the pump head and turn the knurled-head screw 3.1.
2. Rotate the manifold until the inlets of the intake manifold are at the side.
3. Tighten the knurled-head screw.
4. Screw the fastening screws of the eluent line into the inlets of the intake manifold.
5. Use filler caps 3.2 to close unused inlets.

Electrical Connection with the Smartline Pump 1050

For power supply and electric control of the valves of the low pressure gradient module the Smartline Manager 5050 and the Smartline Pump 1050 are connected with a 12-pin connection cable.

In order to correctly install the electrical connection, carry out the following directions in the order given below:

1. Set the Smartline Manager 5050 over the pump.
2. Plug the included connection cable (accessories, Pos. {2.3} in Figure 2) into the "plug for connection" outlet on the rear panel of the pump on the Smartline Manager 5050, Pos. {2.1} in Figure 2.
3. The cable on the rear panel of the Manager 5050 is plugged into the "plug for the PUMP", Pos. {2.1} in Figure 2.

When the Smartline 1050 is turned on, yellow LED's A to D will light up for approximately three seconds and then they will shut off.

A separate 6-/12-pin connection cable is included for the connection of the Smartline Manager 5050 with the Smartline Pump 1000.

Electrostatic discharge can damage the electronics of the pump! Never touch the electric contacts of the Events and Remote terminal strip!

Eluent Connection with the Smartline Pump 1050

The required tube connections for the four LPG connections from the degasser output to the LPG block and from the LPG block to the Smartline Pump pre-cut to fit and included in the accessories (see spare parts and accessories on page 41). Included are 4 PTFE tubes (0.11 m), see Pos. {1.1} in Figure 1, to connect the degasser to the LPG block and 4 PTFE tubes (0.27 m blue marked (D), 0.28 m yellow marked (C), 0.29 m black (B) marked and 0.30 m red marked (A) to connect the LPG block to the pump. The blue PEEK bushings and Tefzel seals must be attached to all of the tubes (see the detailed description in the degasser section of this manual, installation of tube connections on page 24).

Using the 0.11 m tubes connect the degasser to the input of the LPG block "IN" (Pos. {1.1} in Figure 1). Then connect the LPG block "OUT" to the low-pressure mixing block on the pump head (see Pos. {1.6} in Figure 1) using the 0.27 to 0.30 m tubes (see Pos. {1.2} to {1.5} in Figure 1). Close the inputs that are not used with blind fittings.

Please make sure that the solvent bushings and the blind fittings are tight. If they are not, the specifications of the LPG cannot be guaranteed.

Set the desired low-pressure gradient in the pump's setup menu according to SOP 1 and SOP 2.

SOP 1 Selection of the Gradient Mode

This SOP applies to the Smartline Pump 1050 Firmware Revision 1.0X. Use this SOP for selecting a gradient mode.

1. Use the ◀ key to go to the first line.
2. Press the ▶ key approx. 2 seconds to go to the SETUP menu.
3. Press ▲ or ▼ to choose „GRADIENT MODE“.
4. Use the ▶ key to go to the first line.

5. Use the ▲ or ▼ key to choose the desired operation mode. (HPG/LPG).
6. If you select either high or low pressure gradient (i. e. HPG or LPG) you need to specify the channels of used solvents:
LPG: Choose valves A to D
HPG: Choose pump names HPG A to HPG D
7. Move to the second line by pressing ➤ and choose valves A to D (for LPG mode) or pumps actually present (for HPG mode) on the corresponding positions by pressing ▲ or ▼. This selection will apply for programs and define the number and names of solvents used in any gradient.

SOP 2 Selection of operation mode:

This SOP applies to the Smartline Pump 1050 Firmware Revision 1.0X.

1. Press ▲ or ▼ to choose „GRADIENT MODE“.
2. Use the ➤ key to go to the first line.
3. Browse through the GRADIENT MODE options by using the ▲ or ▼ key until LPG appears.
4. Use ➤ to position cursor in the second line.
5. With the ▲ or ▼ keys you can choose ON or -- (off) at any position.
6. Select ON or -- (off) at positions A, B, C and D for the solvent channels in either mode you wish to use. The number of ONs selected refers to all programs in the pumps memory.
7. Press the ➤ key to place the cursor on the next position.
8. By pressing the ➤ key the cursor moves back onto the diamond icon. Now the display should look similar to the graphics in Figure 4.
9. Press the ◀ key approx. 2 seconds to return to the main menu.


Figure 4 Example for a gradient setup with solvents A and B in all programs

The LPG block and the connection tubes should be rinsed with the solvent and filled before they are used.

When the pump is turned off, the LPG block's valves are closed. To rinse the LPG block, set the flow rate to 0 ml/min and set the pertinent mixture portion of the solvent for the respective valve to 100%. After this is completed, switch on the pump, so that the valve can be opened.

Rinse the LPG block with solvent by establishing a connection on the pressure transducer vent capillary on the Smartline Pump 1050 (see pos. {1.3} in Figure 1) using a syringe. Aspirate until no air bubbles appear in the liquid. Repeat this process for all valves to be used.

During operation the yellow "LED's A to D" in front of the Smartline Manager display the respective valve openings.

Vacuum Degasser Unit

Introduction

Working Principle

The Vacuum Degasser is a high-efficiency in-line module that removes dissolved gasses from HPLC solvents. Its unique design assures reliable continuous operation and the highest level of continuous performance available. In the analytical version, up to four solvent lines may be degassed simultaneously by one unit. The low internal volume (285 μl for the analytical version, 7.7 ml for the semipreparative version) of each Teflon AF[®] channel provides for very quick equilibration and very short startup times.

Inside the unit, the solvent flows through a short piece of Teflon AF[®] tubing which is located in a vacuum chamber. Within this chamber a partial vacuum is maintained by a constantly running, low RPM vacuum pump. Dissolved gasses migrate across the tubing wall under a concentration gradient produced by the vacuum as the solvent flows within the coil (see Figure 5). Removed gasses are expelled, and the vacuum is maintained at a constant, preset level by varying the vacuum pump speed as needed.


Figure 5 Principle of gas permeability using a Teflon AF[®] membrane

A special port in the vacuum pump continually flushes the pump head with a small "bleed" of air to remove any solvent vapors which may enter the vacuum pump from the vacuum chamber. This air bleed eliminates the need for any solenoid valves within the system. This patented design results in zero vacuum "hysteresis". Previous designs allowed the vacuum chamber pressure to fluctuate, with the pump cycling on and then off in response to the vacuum level.

Principles of degassing using Teflon AF[®] membranes

This relatively recent addition to the field of degassing has properties not found in other fluoropolymers. The fully amorphous nature of this fluoropolymer and its molecular structure creates a molecular level porosity unlike the mechanically induced porosity in PTFE extruded tubing. In addition, unlike the process used in extruding PTFE, no extrusion agents are needed (like kerosene, etc.) which contaminate mobile phases until they are extracted by the mobile phase over time. Likewise, this molecular structure, combined with the very small surface areas required to degas the mobile phase, reduces the possibility of carryover from one solvent or mobile phase to another to virtually zero.

Teflon AF[®] is so non-polar that it is both solvophobic and hydrophobic. This feature of Teflon AF[®] reduces the possibility of cross-channel contamination from one channel to another, and when combined with the ultra-low internal volumes of Teflon AF[®] channels needed for HPLC flow rates, all but eliminates this cross contamination concern by the

chromatographer. Teflon AF[®] has been used in certain optical systems associated with HPLC for a few years without concern for normal HPLC solvents. However, Teflon AF[®] is soluble in certain solvents (see page 23 Solvents/Mobile Phase) and must not be used to degas these types of solvents.

Teflon AF[®] is permeable to some degree to water vapor whereas PTFE is not. While the vacuum pump in the Vacuum Degasser contains internal provisions for sweeping water or solvent vapor from the pump continuously, it is possible that over time, high concentration buffers may form crystals within the channel due to the loss of water within the channel. The same precautions should be taken to prevent crystallization within these channels as are taken for the HPLC pump. See the "Short-term Shutdown" procedures.

General Description

The Degasser, analytical version, is available for flow rates up to 10 ml/min per channel:

- - Smartline Manager 5050 4 channel analytical degasser unit.

The semipreparative version of the Degasser can be used for flow rates up to 50 ml/min:

- - Smartline Manager 5050 2 channel semipreparative degasser unit.

Connection

Front Panel Connections

Depending upon the model, there are 2 or 4 degassing channels. Pairs of female ¼-28 connectors are located on the front of the Vacuum Degasser module on the Smartline Manager 5050. These are the input and output ports for running up to 4 solvent lines through the Vacuum Degasser. Each channel has an input port and an output port on the same level. In Figure 6 and Figure 7 the connection of both degasser units with LPG unit is shown.

Plugs are provided to seal the ports of unused channels.


Figure 6 *Degasser Tubing Connection, analytical version*


Figure 7 Degasser Tubing Connection, semiprep. Version

Front Panel LED's

The green power LED in the upper left part of the Manager front door indicates the general state of readiness. Two additional LED's are located on the front door of the Manager unit to show the system status of the degasser unit.

POWER (Green)

Indicates when power is applied to the Vacuum Degasser unit (supplied by the Smartline Pump 1050).

Figure 8 Front view of Smartline Manager with Degasser LED's

INIT (Orange)

Indicates when the vacuum level is outside of the acceptable operating range. Normally it will come on at initial power-up and remain lit during pump down. It will go off in a few seconds when the vacuum level goes below 100 mm Hg (133 mbar) absolute. If an error occurs, this LED will flash in one of two modes:

- Flashing on and off in even 1-second intervals: pump was not able to reach vacuum set point, indicating a possible leak in the system.
- Flashing on for 1 second and off for 2 seconds indicates a vacuum signal error. Please contact your Service.

READY (Green)

Indicates when vacuum level is within acceptable operating range. Normally it will come on after the initial pump-down, and remain on as long as the Vacuum Degasser is powered up and vacuum level is below 100 mm Hg absolute.

Rear Panel Connections

Figure 9 Rear view of Smartline Manager with Vacuum Degasser Unit

Exhaust Port

The gas pumped out of the vacuum chamber leaves the unit through the exhaust port (see Figure 9). You can connect a gas tubing to the exhaust port to avoid contaminating the laboratory with solvent gas.

Pump

The degasser unit will be automatically switched on when the Smartline Pump 1050 is switched on. The Smartline Manager must be connected with the pump by a special control cable (G0649).

System Requirements

Solvents/Mobile Phase

Use only HPLC grade solvents. Teflon AF[®] membranes are low in extractables and are very inert, so the degassing System can be used with most solvents including THF and Methylene Chloride. After using saline buffer all lines must be rinsed with water. Do not let the degasser stay for a long period of time with water or buffered soluble solutions. It can lead to a microbacterial growth and/or crystallization. This can result in degasser function defect. Rinse the degasser with isopropanol (always pay attention to the miscibility of the sequential solvents used) if it is not going to be operated for a long period of time (see also page 28 Shut down).

The degassing membrane in the Vacuum Degasser is manufactured from Teflon AF[®]. As with older membranes manufactured from PTFE, Teflon AF[®] is inert to all solvents normally used in HPLC. However, Teflon AF[®] is soluble in perfluorinated solvents such as Fluorinert[®] FC-75 and FC-40 and Fomblin perfluoro polyether solvents from Ausimont. In addition, Freon[®] solvents will adversely affect Teflon AF[®]. Use of such solvents in the Vacuum Degasser will result in the dissolution and hence destruction of the membrane.

Use proper care when handling flammable solvents. Make sure that there are no leaks in the solvent lines. Ensure that hazardous exhaust gases are properly vented.

Chemical compatibility

All parts that contact the mobile phase are made of PEEK, Kel-F[®], Tefzel[®] or Teflon AF[®]. PEEK is not recommended for nitric acid, sulfuric acid and halogenated acid. Organic ethers (THF, MTBE, DMSO) tend to make it swell.

Space Requirements

The Vacuum Degasser Unit is integrated in the Smartline Manager 5050. The Smartline Manager 5050 is designed to sit on a bench top, and is situated in the Smartline LC system between the Smartline eluent tray and Smartline Pump 1050. The Smartline Manager Unit needs a space of 226 x 135 x 390 mm (W x H x D).

Electrical Power Requirements

The AC adapter supplied with the Smartline Pump 1050 incorporates a universal AC input. This allows the instruments to operate at any AC line voltage from 100 to 240 V with a line frequency range of 47 to 63 Hz. With the special gray connection cable (G0649) a regulated voltage of 24V is supplied to the degasser unit.

Installation

Installing the Degasser Unit

The Smartline Manager 5050 with degasser unit can only be used in combination with the Smartline Pump 1050. A special gray connection cable (G0649) must be installed on the back side to ensure the functionality of the degasser unit. To guarantee the electrical power supply of the Smartline Manager 5050 by the Smartline Pump 1050 the devices must be placed in a tower, where the Smartline Manager 5050 is arranged on the Pump 1050 (see Figure 10). As soon as the Smartline Pump 1050 is switched on the vacuum pump of the degasser unit will operate.


Figure 10 Smartline Tower Arrangement

Connecting the Tubings

Solvent lines to be degassed are connected to the Vacuum Degasser's front panel ports, as detailed below. Unused ports must be plugged with the blind fittings to enable the degasser to operate at its peak level of performance.

To make a proper tubing connection, please follow the instructions and scheme. The Smartline Manager 5050 accessory kit is equipped with pre-cut PTFE tubings (see also Table 6). A color identification will help to install the compatible tubings.

All supplied tubings are pre-cut and should be installed to bring about the shortest distance between the connecting ports.

1. Run a blue marked 1 m solvent line of 1/8" O.D. x 1/16" I.D. PTFE chromatography tubing from the solvent supply to the Vacuum Degasser.
2. Install at one end of the solvent tubing the eluent filter.
3. Push the tubing through a PEEK 1/8" male 1/4-28 fitting and slide a ferrule over the tubing end (see Figure 11). Cut the Teflon tubing so the end is flat.

Screw the 1/4-28 fitting into one port on the front of the Vacuum Degasser (channel A with blue marked tubing, for example). Please follow the flow direction on the Degasser unit.
4. The analytical version is to be connected horizontally → and the semipreparative version is to be connected vertically ↑ (comparing to Figure 6 and Figure 7).
5. PEEK connectors should be tightened by hand. Overtightening them will damage the threads.
6. Repeat steps 1 through 3 to connect additional lines to be degassed (channel B with yellow marked tubing, channel C with black marked tubing and channel D with red marked tubing).
7. Once all desired solvent lines have been connected to the Vacuum Degasser, any and all unused ports should be plugged. Use the plugs supplied.
8. Prime each degassing membrane by pulling the solvent from the reservoir through the degassing system. This can be done by connecting a syringe to the tubing or LC pump priming port and drawing air and/or mobile phase into the syringe until no air remains in the tubing, approximately 2 milliliters in case of analytical version and 8 – 10 milliliter in case of semipreparative version.

DO NOT prime the membranes by pushing solvent through the degassing systems. This technique can generate high pressure which might rupture the membrane, even though the Teflon AF[®] membrane is quite rugged. The maximum recommended pressure on the membrane is 1 mPa (100 psig, 7 bar).


Figure 11 Configuration of 1/4-28 Nut, Ferrule and Tubing

Connecting the Vacuum Degasser in a Smartline LPG System

The following picture shows the tubing connections that are typically made between the Vacuum Degasser Unit and other instruments in a Smartline LPG system.


Figure 12 Smartline LPG System Tubing Connections

Operation

Powering up the Vacuum Degasser

1. With the Vacuum Degasser installed into the Smartline System, switch on the Smartline Pump 1050. The green power LED should illuminate. Upon turning on the instrument, the microprocessor immediately examines the vacuum sensor signal to confirm that it is within an expected range. Following the startup test, the microprocessor ramps the vacuum pump to high RPM, to quickly purge the atmosphere from the vacuum chamber. As the vacuum level approaches the preset control value, the pump RPM will slowly ramp down to a low speed (typically 40 to 60 RPM) and will vary slightly as needed under the changing degassing load to maintain a virtually constant vacuum level.

During initial pump-down, the orange INIT LED will be lit. This process only takes 30 seconds. Once the vacuum has reached normal operating level, the orange LED will extinguish and the green READY LED will illuminate.

2. Start solvent flow through the system and check for leaks around the blue flangeless fittings.
If a leak occurs at the connection, tighten the fitting an additional 1/8 turn. If the leak persists, disconnect the leaking fitting and inspect it. If

the nut and ferrule appears to be in good condition, reconnect the fitting. If the leak persists, replace the nut and ferrule and repeat the procedure until you achieve leak-free operation.

The Vacuum Degasser maintains a constant vacuum pressure of 50 mm Hg (67 mbar) absolute (nominal) by varying the speed of the vacuum pump as needed depending on the degassing load in the system.

The pump is designed to run continuously for at least 5 years and is equipped with integral in-pump venting, eliminating the need for stop/start running (U.S. Patent 6,248,157).

The vacuum level and pump speed is constantly monitored by the microprocessor for changes in operating conditions which might be attributed to internal chamber leaks.

If a potential leak is detected (if the pump RPM reaches > 300 for 30 minutes), the pump will be shut down and the orange INIT LED will flash.

The vacuum is maintained as long as the Vacuum Degasser is powered on. Solvent flowing through the Vacuum Degasser will continue to be degassed as long as the instrument is on and running.

3. The Vacuum Degasser Unit will be automatically turned off when the Smartline Pump 1050 to which it is connected is turned off. The vacuum chamber(s) will slowly return to atmospheric pressure when the unit is powered off. This is accomplished by a small, in-line vacuum bleed and reduces the possibility of solvent vapors condensing in the vacuum tubing or pump head.
4. When flushing a solvent line, the single lumen coil inside the chamber contains a very small amount of solvent (approximately 285 microliters for analytical version). When changing from one solvent to another where the final solvent is immiscible with the first, use an intermediate solvent miscible with both the initial and final solvent. Carryover from solvent to solvent is much less than previous PTFE designs. Once air bubbles have been cleared from the solvent line, any further bubbles observed will be coming from the solvent reservoir or from a leaking fitting.

Since there is virtually no solvent retained within the Vacuum Degasser (~285 microliters per channel), priming the system is relatively simple. Using a flow rate of 2 mL/min with the Smartline Pump 1050 for 1-2 minutes ensures that the line from the analytical degasser channel being primed through the proportioning valve on the pump has freshly degassed solvent. If you are using the semipreparative Version the values will be accordingly higher. This dynamic priming method will allow an immediate startup of the analysis upon column equilibration. Contrary to previous PTFE-based degassers, the new Vacuum Degasser, which uses Teflon AF[®] membranes, fully degasses solvents within the time it takes for the volume to pass through the chamber, and yet degasses the solvents as thoroughly as, or better than, PTFE channels containing 40 times more solvent.

Extending the degassing flow rate range

Certain organic solvents used in reversed phase chromatography outgas upon mixing with water, if not properly degassed. These solvents are generally alcohols (e.g. methanol), acetonitrile and tetrahydrofuran. Passing water and methanol through a single channel is generally sufficient to degas these solvents so outgassing does not occur upon mixing when a 75% methanol: 25% water mixture is generated by your

Vacuum Degasser or pump at a flow rate of 5 mL/min. If outgassing does occur, or if a high flow rate is required, it is a general rule that only the organic portion of the mobile phase needs to be passed through a second degassing channel to ensure outgassing does not occur. This is due to the ability of all organic solvents (e.g. methanol) to hold at least 10 times more dissolved atmosphere than water can.

To more thoroughly degas a mobile phase, connect the outlet of the organic channel to the inlet of a second channel and the outlet of the second channel to the pump. This places the two channels in series and doubles the degassing capacity for the organic portion of the mobile phase.

Shutdown

There are two types of shutdown procedures: long-term and short-term.

Short-term Shutdown (Overnight and Weekends)

Observe all precautions pertaining to hazardous solvents and/or those solvents that form harmful deposits or by-products.

1. Flush the column according to the instructions supplied with the column. Flush buffer salts from the system with water. Evaporation leaves salt crystals that may form harmful deposits. Remove chloroform or solvents that can decompose to form hydrochloric acid from the system.

Damage caused by precipitating buffer salts in capillary tubing, or damage resulting from this condition, is specifically excluded from warranty.

2. Remove harmful mobile phases from the Vacuum Degasser and other instruments in the system by carefully purging with pump.
3. After removing harmful mobile phases, prepare the degasser for most mobile phases by flushing it with isopropanol. **To avoid contaminating the system, refill or discard solvents (including water) that were exposed to the environment for more than 24 hours before use.**
4. For weekend storage we recommend flushing 60/40% MeOH/Water through the Vacuum Degasser, pump, column, flow cell (provided your column is compatible with MeOH/Water). Then turn off the Vacuum Degasser, pump and detector.

Long-term Shutdown

1. Follow Short-term Shutdown procedure Steps 1 and 2.
2. Remove the column and direct the pump output tubing to a beaker. Flush the Vacuum Degasser, first with water and then with isopropanol.
3. Turn off the Vacuum Degasser. Then disconnect the tubing between the Vacuum Degasser and solvent reservoirs, and the Vacuum Degasser and pump. Plug all of the ports on the Vacuum Degasser.
4. Store the Vacuum Degasser in a clean, dry location.
5. Before using the Vacuum Degasser, completely purge it with the correct solvent for the column before reconnecting the column and restarting the system.

Maintenance

You can experience long and trouble-free performance from your Vacuum Degasser by performing both routine and preventive maintenance procedures.

Preventative Maintenance

Perform preventive maintenance to ensure that your Vacuum Degasser will perform consistently at an optimal level. To maintain the Vacuum Degasser in the best condition, the following measures are recommended:

1. Adhere to standard laboratory cleanliness practices.
 - Use only high-purity solvents (preferably HPLC grade) for mobile phases. Water should be bottled HPLC grade, or filtered and deionized tap water.
 - Filter all solvents to prevent particulate contamination and tubing blockages.
2. Use only high-purity gases when drying contact areas.
 - Ensure that all new tubing (stainless steel) is passivated and thoroughly flushed before making pump connections.
 - Follow the short- and long-term shutdown procedures that are described above.
3. Routine cleaning of the external surfaces of the instrument can be made using a clean, damp cloth. Immediately clean any spills which occur on or near the instrument using methods appropriate for the type of spill. Some solvents can damage the appearance and function of the instrument.

Routine Maintenance

Routine maintenance is defined as replacing the normal wear items when you notice degradation in performance.

If you have a problem, you can save time and money by referring to the Troubleshooting Guide before calling your Service Representative. Your problem may be minor, and you may be able to correct it yourself using the Troubleshooting Guide to pinpoint the cause.

Never remove the Degasser's cover. There is nothing inside that requires service or maintenance by the customer.

Smart Leak Detection

An additional benefit of maintaining a constant vacuum level is that a potential leak in the vacuum degassing system can be observed by monitoring the RPM of the pump. This "smart leak detection" is a benefit of the patented design of the degasser. If a leak occurs within the chamber, the microprocessor will increase the pump RPM in an attempt to maintain the vacuum level. If the pump cannot maintain the vacuum level (if it runs at an elevated RPM for more than 30 minutes), the yellow LED will flash, indicating a possible leak condition, and the system will shut down and go into a "safe" mode.

Troubleshooting

Table 3 Troubleshooting

Problem	Probable Cause	Solution
Power switch is on, but all 3 LEDs are off, indicating no power to the degasser.	The connection cable (G0649) is not plugged into the pump and manager connector.	Plug the connector cable into the connector.
	Blown fuse.	Contact your technical support.
Orange INIT LED is on steadily, pump is running and RPM seems high.	Pump is in initial pump-down phase or system's degassing demand has increased.	Typically normal operation, although if pump speed continues to rise for an extended period of time (as heard by the pitch of the stepper motor) it could indicate a potential fault condition.
Orange INIT LED is flashing approximately 1 second off, 1 second on. Vacuum pump is not running.	Possible system leak.	Contact your technical support.
Orange INIT LED is flashing approximately 2 seconds off, 1 second on. Vacuum pump is not running.	Possible sensor or control board fault.	Contact your technical support.
Is there a way to check whether the system is operating correctly when Power and Vacuum green LEDs are illuminated, but pump can't be heard running?	Due to the design of the pump and degasser, the pump is virtually silent at low RPM, even though vacuum is good and degassing is normal.	Monitor the UV absorbance of non-degassed methanol at 215 nm versus degassed methanol coming through the degasser. Proper performance of the degasser should decrease the UV absorbance of the methanol significantly.
Bubbles appear through the output tubing.	Loose fitting(s).	Tighten the input and output fittings.
No solvent flow.	Air in the HPLC pump head.	Prime/purge the pump head.
	If a buffer solvent was left in the degasser for some time after use, it may clog the degasser elements.	Try using water to carefully dissolve and rinse out any buffer from the blocked degasser chamber. If this does not help, please contact the technical support department of the manufacturer.

Interface Modul

The Smartline Manager 5050's Interface Module enables the recording or control of up to four instruments that are not supported by the KNAUER software. It can be used for detector recording over the analog inputs, it can control the flow rate for HPLC pumps over the analog outputs, or it can switch valves over the digital outputs. However, the respective instrument must support this option. Please refer to the respective instrument manual regarding this option.

The interface module does not have a stand-alone function; therefore, control software (KNAUER ChromGate® as of version 3.1 or KNAUER ClarityChrom® for Windows) is mandatory.

The interface module functions that are supported by the control software can be found in the software manual. Configuring and programming the Smartline Manager 5050 Interface Module's control software is described in the respective manual of the software.

Instrument Display and Connections

Front view of the instrument


Figure 13 View of the Smartline Manager 5050

Both of the LED's below "Interface" give information about the condition of the interface module. The green "Busy" LED lights shortly during calibration or when a command from the software is received. The instrument is calibrated when it's turned on or a chromatographic run is started. If the control software regularly inquires about the status of the instrument (e.g. in second intervals), then the "Busy" LED will blink permanently.

The red "Error" LED lights when a calibration error is detected during self-test or when the software sends a command to the instrument that does not exist.

Back panel of the instrument


Figure 14 Back panel view of the Smartline Manager 5050 with the Interface Module

Please pay attention to the maximal voltage and current of the inputs and outputs on the interface module and on the instruments that will be controlled. If the maximum voltage and current are exceeded, KNAUER is not liable for defective instruments!

Analog Inputs (AN.IN)

To record from analog signals, up to four channels or instruments can be connected with the interface module. The channels are independent of each other, meaning all channels can be simultaneously operated from either one system or various systems.

Voltage Range:	-2.56 to +2.56 V (bipolar mode)
Absolute maximum ratings:	-10 V to +10 V
Input impedance:	10 M Ω
Maximum resolution:	24 bit
Minimum noise level:	7 μ V (1 Hz, time constant 0.1 s)
Maximum data rate	up to 10 Hz (each channel)
Gain factor (for all channels):	1, 2, 4, 8, 16

Analog Outputs (AN.OUT)

The outputs distribute analog voltage signals, so that instruments containing a respective input can be controlled. Depending on the number of outputs, up to four instruments can be controlled. The outputs can be independently programmed which enables the programming of high pressure gradient system pumps.

Voltage range:	0 to +10 V
Minimum voltage step: (DAC resolution)	2.5 mV (12 bit)
Maximum load resistance:	2 k Ω

Digital Inputs (DIG.IN)

These four inputs are intended as a trigger (start signal) input for the start signal that is distributed by the injection system. They work with contact closure on mass and are compatible with TTL/CMOS or the Open Collector. If the channels are used in separated systems, then every channel can be separately started.

Maximum high level input voltage:	+15 V
Minimum high level input voltage:	+2.5 V
Maximum low level input voltage:	+1 V
Minimum low level input voltage:	-15 V
Maximum input current at V (in) = 0.5 V:	10 mA

Digital Outputs (DIG.OUT)

These four outputs are carried out in electromechanical relays with a single pole double trough (SPDT). They can be selected as NO (normally open; connection 2-3, see the back panel view) or NC (normally closed connection 1-2, see the back panel view). Pulse or continuous signals can be programmed. The length of a pulse contains approximately 1s; the length of the continuous signal can be freely programmed.

Output type:	selectable NO or NC
Maximum switching voltage:	maximum 175 V DC
Maximum switching current:	maximum 0.25 A

RS-232 Interface

This interface is a 9-pin standard RS-232 interface (also known as a serial interface or COM port) and is connected to a PC in order to realize the software control of the interface.

Plug Charging

Pin 2	RxD
Pin 3	TxD
Pin 5	GND

Interface connection settings

Transfer rate:	19200 baud fix,
Data width:	8 bit
Stop bit:	1
Parity:	no

Connecting the Smartline Manager 5050 Interface Module

Power supply

The power supply is secured with the connection of the Smartline Pump 1050 (with a "PUMP" labeled socket on the back panel of the Smartline Manager 5050). The required cable is included in the Smartline Manager 5050's accessory kit.

Smartline Manager 5050 Interface Module on the PC (RS-232)

In order to connect the interface module with a PC, a standard zero modem cable is required and is included in the Smartline Manager 5050 Interface Module's accessory kit. The configuration of the interface are specified under "RS-232 Interface".

Analog Inputs

In order to connect the analog output of the respective instrument (e.g. the integrator output of a detector) with the analog input of the interface module, use the Cinch cable that is made by the instrument's manufacturer and delivered with the instrument. Alternatively, a suitable cable can be found in the Smartline Manager 5050 Interface Module's accessory kit. Please pay attention to the maximal voltage. If it is exceeded, the instrument can be damaged.

The analog input and output of one channel can not be used simultaneously.

Analog Outputs

The analog outputs must be connected with the respective control input of the instrument to be controlled. Information regarding the control input of the instrument can be found in the respective instrument's manual. Please pay attention to the maximal voltage and current of the in and outputs on the interface module and on the instruments that will be controlled. If the maximum voltage and current are exceeded, KNAUER is not liable for defective instruments!

The WAGO plug strips and cable that are required for connection are included in the Smartline Manager 5050 Interface Module's accessory kit.

The analog input and output of one channel cannot be used simultaneously.

Digital Connections

The digital inputs can be used to forward a trigger (start) from the injection unit. The trigger (start) signal must always be connected to the digital input of the same channel to which the detector for data recording is connected.

The digital outputs can be programmed over the control software.

The WAGO plug strips and cable that are required for connection are included in the Smartline Manager 5050 Interface Module's accessory kit.

If a NO is required, pins 2 and 3 (see back panel view) must be connected. However, if a NC is required, pins 1 and 2 must be connected (see back panel view).

If a channel's analog input is used, the digital output (as well as the analog output) of this channel cannot be programmed. In this case, a pulse signal is delivered over the digital output DIG.OUT when a run is started. This can be used as an autozero signal for a detector.

Installation of the WAGO connection cable

In order to clamp the cable on the desired positions of the WAGO plug, please follow the directions below:


Figure 15 Mounting Plug Strips

Example of the Cabling

The exemplary cabling plan demonstrates how the respective connections can be switched. However, it is important that the analog input of a channel is never simultaneously operated with the analog and/or digital outputs of the same channel.


Figure 16 Examples of Wiring

Operation Qualification (OQ)

Only in connection with the pump a detailed function test will be performed. These tests are part of the OQ of the pump.

Installation Qualification (IQ)

Certification on the functionality of the device. During installation of the device, an installation record (IQ document) is created upon request in coordination with the technical service of the manufacturer. This installation protocol needs to be completed in full and be signed by both parties. It serves as proof of the properly executed installation and the functionality of the device.

Maintenance and care

Proper maintenance of your HPLC device will ensure successful analyses and reproducible results.

Technical Support

If you have any technical questions regarding the hardware or software of the manufacturer, please use one of the contact options below:

Hotline of the technical support department:

Languages: German und English
telefonisch erreichbar: 8 – 17Uhr (CET)
Phone: +49-(0)30-809727-0
Telefax: +49(0)30-8015010
E-Mail: info@knauer.net (manufacturer)

Maintenance contract

The following maintenance work on the device may only be performed by the manufacturer or a company authorized by the manufacturer and is covered by a separate maintenance contract:

- Opening the device or removing housing parts.

What maintenance tasks can users perform on the device?

Users may perform the following maintenance tasks themselves:

- Exchange of outer parts subject to wear and tear

Cleaning and caring for the device

Risk of electrical shock or short circuit if cleaning solution enters the device's interior! Only lightly moisten the cleaning cloth!

All smooth surfaces of the device can be cleaned with a mild, commercially available cleaning solution, or with isopropanol.

Environmental protection

Disposal

Drop off old devices at the local municipal waste facilities or send back old devices to the manufacturer where it will be disposed of properly.

Decontamination

Contamination of devices with toxic, infectious or radio-active substances poses a hazard for all persons during operation, repair, sale and disposal of a device.

Danger caused by toxic, infectious, or radio-active substances! A contaminated device must never be submitted for repairs, sold, or disposed of! Contract a specialist company to decontaminate the device or perform the decontamination yourself if you have the required expertise!

All contaminated devices must be properly decontaminated by a specialist company or the operating company before they can be recommissioned, repaired, sold, or disposed of.

All materials or fluids used for decontamination must be collected separately and disposed of properly.

Storage

Ambient storage conditions for the device

Temperature range: 4–40 °C; 39.2–104 °F

Humidity: Below 90 % humidity (non-condensing)

Technical data

Table 4 Technical Data and Specifications

Specification	Description
Dimensions	226 x 135 x 390 mm (W x H x D)
Weight	Max. 4.8k g (complete Smartline Manager 5050)
Channels	1-4 independent / 1-2 independent
Degassing Flow Path I.D.	1.14 mm analytical version / 1.77 mm semipreparative version
Degassing Process	Gas permeation through a fluoropolymer membrane
Maximum Flow Rate	10 ml/min analytical / 50 ml/min semipreparative
Degassing Capacity	~2 ppm at 1 ml/min for analytical version
Degassing Efficiency	< 0.5 ppm residual dissolved O ₂ at 1 ml/min
Degassing Channel Pressure Rating	70 PSIG
Dead Volume	~285 µl per channel for analytical version ~7.7 ml per channel for semipreparative version
Optimized Flow Rate (for isocratic or gradient formed 50:50 methanol/water)	3 ml/min analytical version / 30 ml/min preparative version
Material in contact with solvents	PEEK, Glass-filled PTFE, Teflon AF [®]
Solvent Applicability	universal, except hydrochloric acid and halogenated carbohydrates especially hexafluoro isopropanol (HFIP)
Power	via Smartline Pump 1050
Supply voltage range	100-240 V
Supply frequency	50–60 Hz

Specification	Description
Power consumption	Maximum 40 W
IP protection class	IP-20
Operating Temperature	4 to 40 °C
Operating Relative Humidity	20 to 80 % RH (without condensation)
Altitude	0 to 2000 Meters
Storage temperature	-20 to +60 °C
Storage Relative Humidity	20 to 80% RH (without condensation)

Table 5 Technical Data and Specifications of NDG

Specification	Description
Eluent count	Up to 4
Maximum flow per channel	50 ml/min
Gradient range	0 – 100 %
Gradient gradation	1 % steps
Control	Through the Smartline Pump 1050 or ChromGate® or ClarityChrome® software
Material in contact with solvent standard	Stainless steel, FFKM, PEEK, Tefzel
Material in contact with solvent Bio	FFKM, PEEK, Tefzel

Delivery program

The tables listed below contain all necessary positions for the successful installation of the Smartline Manager 5050 unit. The tubings are colour highlighted to make the installation and allocation easier.

Table 6 List of Smartline Manager 5050 Accessory A5331

Description	Order Number
Syringe	N0102
Stainless steel solvent filter, 2x	A1149V1
PTFE tubing 1 m with eluent filter (red – channel D)	G1440-1
PTFE tubing 1 m with eluent filter (black – channel C)	G1440-2
PTFE tubing 1 m with eluent filter (yellow – channel B)	G1440-3
PTFE tubing 1 m with eluent filter (blue – channel A)	G1440-4
PTFE tubing 3.2 mm (1/8") OD	A0732*
Flangeless ferrule 1/8, 3 x 10 pieces	M1057
Flangeless fittings 1/8, 3 x 10 pieces	M0689
Luerlock needle 1.5 x 50 mm	M1551
PTFE tubing 0.27 m (from LPG-channel D to pump) red	G0766
PTFE tubing 0.28 m (from LPG-channel C to pump) black	G0767
PTFE tubing 0.29 m (from LPG-channel B to pump) yellow	G0768
PTFE tubing 0.30 m (from LPG-channel A to pump) blue	G0769
Connecting cable 12 pole	G1708
WAGO 9 pole	G0750
Null Modem cable	M1922
Connector cable with wire ends	G1021
Flat band cable 10-pole	A1467
Complete connection 4x gradient	G4145
PTFE tubing 110 mm	P9170
Smartline Manager 5050 Manual	V7599

* A0732 PTFE tubing, Length 3 m (can be pre-cutted by the user)

Table 7 List of Smartline Manager 5050 Accessory A5332

Description	Order Number
Syringe	N0102
Stainless steel solvent filter, 2x	A1149V1
PTFE tubing 1 m with eluent filter (red – channel D)	G1440-1
PTFE tubing 1 m with eluent filter (black – channel C)	G1440-2
PTFE tubing 1 m with eluent filter (yellow – channel B)	G1440-3
PTFE tubing 1 m with eluent filter (blue – channel A)	G1440-4
PTFE tubing 3.2 mm (1/8") OD	A0732*
Flangeless ferrule 1/8, 3 x 10 pieces	M1057
Flangeless fittings 1/8, 3 x 10 pieces	M0689
Luerlock needle 1.5 x 50 mm	M1551
PTFE HPG tubing 250 mm	G6113
PTFE HPG tubing 450 mm	G6112
Connecting cable 12 pole	G1708
WAGO 9 pole	G0750
Null Modem cable	M1922
Connector cable with wire ends	G1021
Flat band cable 10-pole	A1467
Smartline Manager 5050 Manual	V7599

* A0732 PTFE tubing, Length 3 m (can be pre-cutted by the user)

Table 8 List of Smartline Managers 5050 Accessory A5333

Description	Order Number
Syringe	N0102
Stainless steel solvent filter, 2x	A1149V1
PTFE tubing 1 m with eluent filter (red – channel D)	G1440-1
PTFE tubing 1 m with eluent filter (black – channel C)	G1440-2

Description	Order Number
PTFE tubing 1 m with eluent filter (yellow – channel B)	G1440-3
PTFE tubing 1 m with eluent filter (blue – channel A)	G1440-4
PTFE tubing 3.2 mm (1/8") OD	A0732*
Flangeless ferrule 1/8, 3 x 10 pieces	M1057
Flangeless fittings 1/8, 3 x 10 pieces	M0689
Luerlock needle 1.5 x 50 mm	M1551
PTFE tubing 110 mm	P9170
PTFE tubing 0.27 m (from LPG-channel D to pump) red	G0766
PTFE tubing 0.28 m (from LPG-channel C to pump) black	G0767
PTFE tubing 0.29 m (from LPG-channel B to pump) yellow	G0768
PTFE tubing 0.30 m (from LPG-channel A to pump) blue	G0769
Connecting cable 12 pole	G1708
Complete connection 4x gradient	G4145
Smartline Manager 5050 Manual	V7599

* A0732 PTFE tubing, Length 3 m (can be pre-cutted by the user)

Table 9 List of Smartline Managers 5050 Accessory A5338

Description	Order Number
Syringe	N0102
Stainless steel solvent filter	A1149V2
Luerlock needle 1.5 x 50 mm	M1551
PTFE tubing 3.2 mm 1/8 AD	A0873
PTFE tubing 3.2 mm (1/8") OD	A0732*
Flangeless ferrule 1/8, 3 x 10 pieces	M1057
Flangeless fittings 1/8, 3 x 10 pieces	M0689
Connecting cable 12 pole	G1708
WAGO 9 pole	G0750

Description	Order Number
Null Modem cable	M1922
Connector cable with wire ends	G1021
Flat band cable 10-pole	A1467
Smartline Manager 5050 Manual	V7599

* A0732 PTFE tubing, Length 3 m (can be pre-cutted by the user)

Table 10 List of Smartline Managers 5050 Accessory A5331-1, A5333-1, A5338-1

Description	Order Number
Syringe	N0102
PTFE tubing 1 m with eluent filter (red – channel D)	G1440-1
PTFE tubing 1 m with eluent filter (black – channel C)	G1440-2
PTFE tubing 1 m with eluent filter (yellow – channel B)	G1440-3
PTFE tubing 1 m with eluent filter (blue – channel A)	G1440-4
PTFE tubing 110 mm	P9170
PTFE tubing 140 mm	P9171
PTFE tubing 440 mm	P9172
Flangeless ferrule 1/8, 3 x 10 pieces	M1057
Flangeless fittings 1/8, 3 x 10 pieces	M0689
Luerlock needle 1.5 x 50 mm	M1551
PTFE tubing 0.27 m (from LPG-channel D to pump) red	G0766
PTFE tubing 0.28 m (from LPG-channel C to pump) black	G0767
PTFE tubing 0.29 m (from LPG-channel B to pump) yellow	G0768
PTFE tubing 0.30 m (from LPG-channel A to pump) blue	G0769
PTFE HPG tubing 250 mm	G6113
Connecting cable 12 pole	G1708
Complete Connection 4x gradient	G2783
Intake manifold, 4x way	G2783V1
Pipe connection, titanium	G2779V2

Description	Order Number
Capillary pipe, titanium	P9090V2
Dynaseal screw fitting for PEEK sample loop	P0968
Cutting ring 1/16"	P0307
Split-grooved clamping ring 1/16"	P0308
Sealing ring for 1/16" PEEK	P0506
Connection piece PEEK, 4-way	G2747
Screw fitting, outlet side	P3169
Spacer sleeve	P3349
Cylinder screw ISO 4762- M4x60 A2	R0534
Connection cable DIN plug on DIN plug	G0649
Smartline Manager 5050 Manual	V7599

Table 10 List of Smartline Manager 5050

Description	Order Number
Power connector	M1479
WAGO 9 pole	G0750
Null Modem cable	M1702-1
Connector cable with wire ends	G1021
Flat band cable 10-pole, 1.5	A1467

Table 11 List of Smartline Manager 5050

Description	Order Number
Smartline Manager 5050: Interface module, degasser module analytical, LPG module	A5331
Smartline Manager 5050: Interface module, degasser module analytical, LPG module, biocompatible	A5331-1
Smartline Manager 5050: Interface module, degasser module analytical	A5332
Smartline Manager 5050: Degasser module analytical, LPG module	A5333
Smartline Manager 5050: Degasser module analytical, LPG module, biocompatible	A5333-1
Smartline Manager 5050: Interface module, degasser module semi- preparative, LPG module	A5338

Description	Order Number
Smartline Manager 5050: Interface module, degasser module semi-preparative, LPG module, biocompatible	A5338-1
Smartline Manager 5050 Manual	V7061

Legal information

Warranty conditions

The factory warranty for the device is valid for 12 months after the date of dispatch. All warranty claims shall expire in the event that any unauthorized changes are made to the device.

During the warranty period, any components with material or design-related defects will be replaced or repaired by the manufacturer free of charge.

This warranty excludes the following:

1. Accidental or willful damage
2. Damage or errors caused by third parties that are not contractually related to the manufacturer at the time the damage occurs
3. Wear parts, fuses, glass parts, columns, light sources, cuvettes and other optical components
4. Damage caused by negligence or improper operation of the device and damage caused by clogged capillaries
5. Packaging and transport damage

In the event of device malfunctions, contact:

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The packaging of our devices provides the best possible protection against transport damage. However, immediately inspect each delivery for signs of transport damage. If the shipment is incomplete or damaged, inform the manufacturer within three workdays. Also inform the freight carrier about transport damage.

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Smartline Manager 5050 Order numbers:

A5331, A5331-1, A5332, A5333, A5333-1, A5338, A5338-1

Complies with the following requirements and product -specifications:

- DIN EN 60799 (June 1999) Electrical accessories – Cord sets and interconnection cord sets
- DIN EN 61010-1 (August 2002) Safety requirements for electrical equipment for measurement, control and laboratory use
- Low voltage directive (2006/95/EC)
- DIN EN 61000-3-2 (October 2006) Electromagnetic compatibility (EMC) Part 3-2
- EMC directive (2004/108/EC)
- DIN EN 61326-1 (October 2006) Electrical equipment for measurement, control and laboratory use – EMC requirements
- Directives for an environmentally sound use of electrical and electronic equipment
- RoHS directive 2002/95/EC (February 2003) on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- WEEE directive 2002/96/EC (February 2003) on waste electrical and electronic equipment

Berlin, 2014-03-04

Dr. Alexander Bünz (Managing Director)

The mark of conformity has been applied to the rear panel of the device.

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© Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
All rights reserved.
The information in this document is subject to
change without prior notice. Translation of the
original German edition of this manual, version 1.2
2014-03-04
Printed in Germany on environmentally friendly
paper from sustainable forests.

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