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CONTENTS

Using this Manual	5
SOPs in this manual	5
General	7
Capillary Connections.....	7
Manually driven 6-Port-3-Channel Valves.....	9
Delivery Content	9
Layout and Function	10
Installation.....	11
Operation	12
Using the valve as a Switching Valve.....	12
Electrically driven 6-Port-Multi-Channel Valves	13
Delivery Content	13
Layout and Function	14
Use as an Injection Valve.....	15
Using the valve as a Switching Valve.....	15
Special version for the 6-Port-2-Channel Valve	15
Manual 6-Port-2-Channel Injection Valve.....	16
Delivery Content	16
Assembly and Function.....	16
Manual 6-Port-3-Channel Switching Valves.....	18
Delivery Content	18
Layout and Function	19
7-Port-1-Channel Valves	19
Delivery Content	19
Layout and Function	20
Examples for using the 7-Port-1-Channel Switching Valves	20
Step-function gradient	20
Fraction collector	21
Column switching	21
17(13)-Port-1-Channel Switching Valves	22
Delivery content	22
Layout and Function	22
Examples for switching with Multiport Valves.....	23
Column selection with a 2-Channel Valve	23
Using a precolumn to precondition the sample	24
Selecting a column with a 3-Channel Valve	24
Precolumn back-flushing	24
Column back-flushing	25
Sample enrichment.....	25
Alternating enrichment with a precolumn through rinsing.....	26
Alternating enrichment with precolumn back-flushing	26
Alternating sample preconditioning with precolumn through rinsing.....	27
Alternating sample preconditioning with precolumn back-flushing.....	27
Repair and Service Note	28
Accessories and Spare Parts	30
Spare Parts.....	30
Accessories	30
Guarantee statement.....	58
Declaration of conformity.....	59

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Using this Manual

Special Warnings are indicated by the marginal warning sign and printed in bold letters.

Important Hints are marked by the marginal hand symbol.

The marginal **light bulb** symbol indicates **helpful advice**.

SOPs in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks for operating Injection and Switching Valves. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and can be copied, applied, signed, and filed in order to document the performance of the instrument.

It is very important to follow all instructions and SOPs in this manual when operating the valves. This ensures proper results and longevity of your equipment.

Table 1 List of SOPs in this manual

SOP 1	Capillary connection using the DYNASEAL system	8
SOP 2	Capillary connection using the UNF system	8
SOP 3	Connections to a 6-Port-3-Channel Injection Valve	12
SOP 4	Sample Injection	12
SOP 5	Mounting the valve to the electrical drive	14
SOP 6	Installing the valve into the thermostat	17
SOP 7	Connections on a 6-Port-2-Channel Injection Valve	17
SOP 8	Dismantling the valve body	28
SOP 9	Inserting a new needle seal in the injection port	29

This manual refers to the KNAUER Injection and Switching Valves, listed in the following table

Table 2 List of KNAUER Injection and Switching Valves

Order No.	Ports	Channels	Function	Material	Connection	Drive
A1357	6	3	injection switch	stainless steel	1/16"	manual
A1358	6	3	injection switch	PEEK	1/16"	manual
A1359	6	3	injection switch	stainless steel	1/8"	manual
A1360	6	3	injection switch	PEEK	1/8"	manual
A1361	7	1	switch	stainless steel	1/16"	manual
A1362	7	1	switch	PEEK	1/16"	manual
A1363	7	1	switch	stainless steel	1/8"	manual
A1364	7	1	switch	PEEK	1/8"	manual
A1365	6	2	injection	stainless steel	1/16"	manual ¹⁾
A1366	6	2	injection	PEEK	1/16"	manual ¹⁾
A1367	6	3	switch	stainless steel	1/16"	manual ¹⁾
A1368	6	3	switch	PEEK	1/16"	manual ¹⁾
A1369	6	3	injection switch	stainless steel	1/16"	K-6
A1369-1	6	2	switch	stainless steel	1/16"	K-6
A1370	6	3	injection switch	PEEK	1/16"	K-6
A1371	6	3	injection switch	stainless steel	1/8"	K-6
A1371-1	6	2	injection switch	stainless steel	1/8"	K-6
A1372	6	3	injection switch	PEEK	1/8"	K-6
A1373	7	1	switch	stainless steel	1/16"	K-6
A1374	7	1	switch	PEEK	1/16"	K-6
A1375	7	1	switch	stainless steel	1/8"	K-6
A1376	7	1	switch	PEEK	1/8"	K-6
A1378	13	1	switch	stainless steel	1/8"	K-12
A1379	17	1	switch	stainless steel	1/16"	K-16
A1379-1	17	1	switch	stainless steel	1/16"	K-16

¹⁾ built-in version for thermostat A0585

General

All KNAUER Injection and Switching Valves were developed for HPLC use. The valves are available in stainless steel and PEEK. If the PEEK version is used, the sample and eluent only come into contact with the polymer. Thus, the prerequisite for working under inert conditions is met.

In general, the valves are adjusted for pressure up to 300 bar. Upon request, they can be adjusted for pressure up to 400 bar. It is also possible for our service department to readjust the valves.

The Multi-Channel Injection and Switching Valves are equipped with a reed contact and magnet. By switching the valve to the INJECT position a contact is made. Switching the valve to the LOAD position will reopen the same contact.

The reed contact is mainly used as a start signal for either the HPLC software or an integrator. It is also used with detectors to correct the baseline with each injection by performing an auto zero.

If a pulse is needed instead of a permanent signal, a pulser (Order No. A1137) can be used. It generates a pulse when switching from LOAD to INJECT, but not when switching from INJECT to LOAD.

Further details regarding the reed contact can be found in the sections on the individual valve groups or in the respective manuals for the devices.

Do not use the valves in dry conditions. Otherwise, it can result in a decrease in the lifetime of the rotor sealing. This is especially true of the motor driven valves.

If aqueous salt or buffer solutions are used as eluents, it is strongly recommended to rinse the valve with distilled water and/or methanol after use. This will prevent the formation of crystalline deposits inside the valve.

Make sure the valve is switched quickly to avoid the pump from being blocked by pressure surges!

Capillary Connections

In general, the capillary connections to the Injection and Switching Valves can be made with stainless steel or PEEK capillaries. In both cases the capillaries can be fixed with the DYNASEAL connection system. However, stainless steel capillaries provide a longer secure connection with the enclosed UNF bushings.

To avoid damage do not use steel connections for PEEK valves. Only use polymer capillaries, split-grooved clamping rings, polymer sealing rings, and DYNASEAL bushings. The DYNASEAL bushings should be tightened by hand only.

All capillary connections must be made with minimal dead volume. Therefore, use the shortest capillaries possible with a small inner diameter.

Fig. 1 shows a schematic diagram of the UNF and DYNASEAL connection system.


Fig. 1 Principle of the UNF and DYNASEAL connection systems

Both systems are available with short and long bushings. The shape of the bushings is not relevant to their function. Only make your selection according to the local geometrical conditions.

The accessibility of the bushings can often be optimized by alternately inserting the bushing types.

SOP 1 Capillary connection using the DYNASEAL system

1. Make sure that the capillary is squarely cut off. If necessary, cut the capillary using the tube cutter (Order No. 0569).
2. Push the bushing, the split-grooved clamping ring, and the sealing ring onto the capillary. Pay attention to the sequence and alignment of the fittings.
3. Push the capillary into the valve port as far as possible.
4. Hand tighten the bushing while pushing the capillary towards the port.

SOP 2 Capillary connection using the UNF system

1. Make sure the capillary is squarely cut off. If necessary, cut the capillary using the cutting pliers (Order No. 0809) and deburr it with the capillary cleaning set (Order No. A0137).
2. Push the bushing and the ferrule onto the capillary. Pay attention to the sequence and alignment of the fittings (see Fig. 1).
3. Push the capillary into the valve port as far as possible.
4. Hand tighten the bushing while pushing the capillary towards the port.
5. Slightly tighten the bushing with a hexagonal spanner.


Fig. 2 Seating of capillary connectors

Manually driven 6-Port-3-Channel Valves

A1357, A1358, A1359, A1360

The valves are pressure stable up to 300(400) bar. Pay attention to the general introductory sections!

Delivery Content

Order No.	Pieces	Item	Details
A1357	1	6-Port-3-Channel Injection Valve	Stainless steel, 1/16"
	4	Ferrule	1/16" 1.4401
	2	Double-cone sealing	1/16" PETP
	1	Bushing UNF 10/32	Stainless steel, 12mm
	3	Bushing UNF 10/32	Stainless steel, 15mm
	2	Open-end wrench	1/4" and 5/16"
	2	Needle seal 34mm	AD1.6 x ID 0.9mm
	1	DIN911 Allen-wrench	No.3, nickel plated
	1	Manual	
	A1358	1	6-Port-3-Channel Injection Valve
1		Double-cone sealing	1/16" PETP
4		Split-grooved clamping ring	1/16"
4		Polymer sealing ring	For 1/16" PEEK
4		DYNASEAL bushing	1/16" long
2		Needle seal 34mm	AD1.6xID0.9mm
1		DIN911 Allen-wrench	No.3, nickel plated
1		Manual	
A1359	1	6-Port-3-Channel Injection Valve	Stainless steel, 1/8"
	6	Bushing	1/8" M8x1 SW10
	6	Ferrule	1/8"
	2	Open wrench	8/10"
	1	DIN911 Allen-wrench	No.3, nickel plated
	1	DYNASEAL bushing	1/8" length M8x1
	1	Double-cone sealing	1/8" PETP
	2	Needle seal 1/8" 37mm	AD1/8" x ID1.5mm
	1	Injection syringe	10ml, Luer connector
	1	Luer-Lock cannula	1.5x50, straight
	1	Manual	
A1360	1	6-Port-3-Channel Injection Valve	PEEK, 1/8"
	6	DYNASEAL bushing	1/8" length M8x1
	6	Polymer seling ring	For 1/8" PEEK
	1	Needle seal 1/8" 37mm	AD1/8"xID1.5
	6	Split-grooved clamping ring	1/8"
	1	Double-cone sealing	1/8" PETP
	1	Injection syringe	10ml, Luer connector
	1	Luer-Lock cannula	1.5x50, straight
	1	DIN911 Allen-wrench	No.3, nickel plated
	1	Manual	

Layout and Function

The four, manually driven, 6-Port-3-Channel Valves A1357, A1358, A1359, and A1360 only differ in material (stainless steel and PEEK) and/or in the port, and some internal dimensions. Therefore, they have been described together.

The body of the valves within the housing consists of a rotor above automatically-adjusting spring disks, start up disks, and an axial needle bearing. Connected to the rotor is a lever which moves perpendicular to the rotor axis and is used to move the rotor. The lever cannot be moved through more than 60°. Within the housing there is a reed contact for determining the position of the rotor. This is activated by a magnet located in the lever. The rotor seal is located on the rotor.

Fig. 3 Manually operated 6-Port-3-Channel Valve


Fig. 4 Manually operated 6-Port-3-Channel Valve, exploded view

The rotor can be turned by 60° with the lever.

Lever left:	LOAD position
Lever right:	INJECT position

In the INJECT position the magnet activates the reed contact. Making this contact can start a gradient programmer or detector. Different adjacent ports are connected via the passages in the rotor sealing depending on the position of the rotor (see Fig. 5).


Fig. 5 Port connections in the LOAD and INJECT positions

In the LOAD position the sample loop will be filled. When switched to the INJECT position the sample will be injected. The eluent passes the sample loop in opposite direction to minimize the band broadening. This is especially advantageous if sample loops are partially filled.


Fig. 6 Flow directions while loading (sample) and injecting (eluent)

Installation

The Universal Mounting Bracket (A4018) can be used to fix the valve to any of the KNAUER pumps. Use the long size holes on the short bracket side and the tapholes on the side of the pump to fix the valve with M4 screws.

For the small HPLC pump K-120 an especially suited angel, (A0406), can be used.


Fig. 7 Universal Mounting Bracket for valve attachment

Use M4 screws to secure the valve with the tapholes on the rear panel of the valve (1) Fig. 7. The feedthroughs (2) are used for capillaries and position (3) is where the column holder is secured.

Connecting capillaries to the valve is described in the general section, *Capillary Connections* on page 7. The port allocation depends on how the valve is used. Therefore, it is described in the following section.

Operation

Using the valve as an Injection Valve

To use the valve as an Injection Valve a syringe port consisting of a normal DYNSEAL bushing and a needle seal (see Fig. 32 on page 29) has to be inserted into port 1.

SOP 3 Connecting a 6-Port-3-Channel Injection Valve

1. Insert the **syringe port** into **port 1** and gently tighten the bushing.
2. Insert an injection syringe into the needle seal.
3. Together with the syringe, push the the needle guide to the end of the injection port.
4. Connect a sample loop to **port 2** and **port 5**.
5. Connect the **pump** to **port 4**.
6. Connect **port 6** to the waste container with teflon tubing or a capillary.

Try to keep the overflow at approximately the same height as the syringe port. This will avoid a suction effect caused by the difference in height.

7. Rinse the valve with the eluent.
8. After rinsing, connect the column to port 3.

Injecting a Sample

Only use a square (90°) Luer-Lock injection syringe 1/16" diameter without ridges.

The size of the syringe and the injection volume should correspond.

A minimum of a 5-fold dose of sample solution should be used for the sample injection if working with full loop injection. This results in a relative deviation of $\leq 0.2\%$.

The sample loop can also be partially filled. Then the sample volume is determined by the injected volume.

SOP 4 Sample Injection

1. Set the valve to LOAD.
2. Fill the syringe with the sample and remove the air bubbles.
3. Insert the syringe into the syringe port in port 1 as far as possible.
4. Empty the syringe contents into the loop and leave the syringe in position.
5. Set the valve to INJECT and (if necessary) simultaneously start the data acquisition.
6. Set the valve to LOAD and remove the syringe from the port.
7. If necessary, flush the inlet port and sample loop with eluent (Generally, only required for trace analysis since contamination is negligible).

Make sure the valve is switched quickly to avoid the pump from being blocked by pressure surges!

Using the valve as a Switching Valve

The four manually driven 6-Port-3-Channel Valves can also be used for various switching tasks. Details can be found in the *Examples on Switching with Multiport Valves* section on page 23. In detail descriptions are given for the switching possibilities of 2-Channel and 3-Channel Valves. In this respect, differentiation is not necessary between manually and electrically driven valves.

Electrically driven 6-Port-Multi-Channel Valves

A1369, A1370, A1371, A1372

The valves are pressure stable up to 300(400) bar. Pay attention to the general introductory sections!

Delivery Content

Order No.	Quantity	Article	Details
A1369 A1369-1	1	6-Port-3-Channel Injection Valve (6-Port-2-Channel Valve)	Stainless steel, 1/16"
	8	Ferrule	1/16"
	2	Biconical sealing ring	1/16" PETP
	3	Bushing UNF 10/32	Stainless steel, 12mm
	3	Bushing UNF 10/32	Stainless steel, 15mm
	2	Open-end wrench	1/4" and 5/16"
	2	Needle seal 34mm	AD1.6 x ID 0.9mm
	1	Manual	
A1370	1	6-Port-3-Channel Injection Valve	PEEK, 1/16"
	1	Biconical sealing ring	1/16" PETP
	4	Split-grooved clamping ring	1/16"
	4	Sealing ring	For 1/16" PEEK
	4	DYNASEAL bushing	1/16" long
	2	Needle seal 34mm	AD1.6xID0.9mm
	1	DIN911 6-hexagonal key	No.3, nickel plated
1	Manual		
A1371 A1371-1	1	6-Port-3-Channel Injection Valve (6-Port-2-Channel Valve)	Stainless steel, 1/8"
	6	Bushing	1/8" M8x1 SW10
	6	Ferrule	1/8"
	2	Open-end wrench	8/10"
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	DYNASEAL bushing	1/8" length M8x1
	1	Dual-cone sealing ring	1/8" PETP
	2	Needle seal 1/8" 37mm	AD1/8" x ID1.5mm
	1	Disposable syringe	10ml, Luer connection
	1	Luer-Lock-cannula	1.5x50, dull
1	Manual		
A1372	1	6-Port-3-Channel Injection Valve	PEEK, 1/8"
	6	DYNASEAL bushing	1/8" length M8x1
	6	Sealing ring	For 1/8" PEEK
	1	Needle seal 1/8" 37mm	AD1/8"xID1.5
	6	Split-grooved clamping ring	1/8"
	1	Dual-cone sealing ring	1/8" PETP
	1	Disposable syringe 10ml	Luer connection
	1	Luer-Lock-cannula	1.5x50, dull
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	

Layout and Function

The four electrically operated 6-Port-3-Channel Valves are only different in material (i.e. stainless steel or PEEK) and/or in the port, and a few internal dimensions. Therefore, they have been described in the same category.

The construction and function corresponds with that of the manually operated valve (see page 10). The main difference is that the manual valve does not have an operating lever. Instead the valves are moved by the Smartline Valve Drive 6 (Order No. A55020) or the Electrical Valve Drive K-6 (Order No. A0377).

Fig. 8 7-Port-1-Channel Valve with Valve Drive K-6

SOP 5 Mounting the valve to the electrical drive

1. Check the positions of the DIP switches on the Smartline Valve Drive 6 or Electrical Valve Drive K-6, described in the Valve Drive manual and correct them accordingly.
2. Screw the adaptor plate with both of the M4-flat-head screws on the back panel of the valve
3. Set the drive to position 1 or LOAD.
4. Push the valve on the driving axle located on the front side of the instrument.
5. Make sure the mark on the valve matches the one above the wheel. The valve position, 1, must be on top.
6. Secure the valve and adapter plate with the two mounting screws onto the housing.

All 6-Port Valves can only be switched with an Smartline Valve Drive 6 or Electrical Valve Drive K-6. Operation with a 12 or 16 positions drive results in mechanical malfunctions and may damage the HPLC system.

Fig. 9 illustrates the connection of the capillary connections to the individual devices as well as the eluent flow through the valve.


Fig. 9 Port connection in the LOAD and INJECT position

Use as an Injection Valve

Operating the electrical moving injection valve corresponds to that of the manually operated valve.

Production of the capillary connections can be found under **SOP 3** *Connecting a 6-Port-3-Channel Injection Valve* on page 12.

Injecting a sample

Only use a square (90°) Luer-Lock injection syringe 1/16" diameter without ridges.

The size of the syringe and the injection volume should correspond.

A minimum of a 5-fold dose of sample solution should be used for the sample injection if working with full loop injection. This results in a relative deviation of $\leq 0.2\%$.

The sample loop can also be partially filled. Then the sample volume is determined by the injected volume.

For more information on sample injection see **SOP 4** *Sample Injection* on page 12.


Fig. 10 Feed pump connection without sample loop

Using the valve as a Switching Valve

The four manually driven 6-Port-3-Channel Valves can also be used for various switching tasks. Details can be found in the *Examples on Switching with Multiport Valves* section on page 23. In detail descriptions are given for the switching possibilities of 2-Channel and 3-Channel Valves. In this respect, differentiation is not necessary between manually and electrically driven valves.

Special version for the 6-Port-2-Channel Valve

The 3-Channel Valves (A1369 and A1371) can also be used as a 2-Channel Valves (A1369-1 and A1371-1). In order to do so disassemble the rotor seal and rotate it 180° diagonally across the rotor axle before reinstallation. The delivery content corresponds to that of the 3-Channel Valves.


Fig. 11 Internally connected ports according to the valve position

Manual 6-Port-2-Channel Injection Valve

A1365, A1366

The valves are pressure stable up to 300(400) bar. Pay attention to the prefixed sections related to all valves!

Delivery Content

Order No.	Quantity	Article	Details
A1365	1	6-Port-2-Channel Injection Valve	Stainless steel, 1/16"
	4	Ferrule	1/16"
	4	Bushing UNF 10/32	Stainless steel, 12mm
	2	Open-end wrench	1/4" and 5/16"
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	
A1366	1	6-Port-2-Channel Injection Valve	PEEK, 1/16"
	4	DYNASEAL bushing	1/16" long
	4	Split-grooved clamping ring	1/16"
	4	Sealing ring	For 1/16" PEEK
	1	Manual	

Assembly and Function

Both of the manually operated 6-Port-2-Channel Valves were designed specifically for installation in the columns thermostat (Order No. A0585). They are only differentiated by the material (stainless steel or PEEK) and are, therefore, described together.


Fig. 12 6-Port-2-Channel Injection Valve

The most significant difference in the 6-Port-3-Channel Valves is that the injection port is exported through the elongated rotor out of the thermostat. For this reason a longer injection port must be used (see Fig. 13).


Fig. 13 Exploded view of a 2-Channel Injection Valve

SOP 6 Installing the valve into the thermostat

1. Unscrew the injection port
2. Release the lever from the valve by removing both screws.
3. Push the rotor axle from inside through the notch in the top side of the thermostat.
4. Secure the valve on the top side of the thermostat with two M4°screws.
5. Remount the lever on the axle.
6. Screw the injection port into the hole where the lever axle is located.
7. Push the reed contact connection outwardly through a capillary feedthrough.

SOP 7 Connections on a 6-Port-2-Channel Injection Valve

1. Insert the syringe port into the outer port of the valve and gently tighten the bushing by hand.
2. Connect the sample loop diagonally to port 3 and 6.
3. Connect the pump to port 4.
4. Ports 1 and 2 are connected by teflon tubing or capillaries with a waste container.

Try to keep the overflow at approximately the same height as the syringe port. This will avoid a suction effect caused by the difference in height.

5. Rinse the valve with the eluent.
6. After rinsing, connect the column to port 5.

The connections from the pump to the detector and waste container are guided through the capillary feedthrough located on the side of the thermostat.

Fig. 14 illustrates the connection of the capillary connections to the individual devices such as the eluent flow through the valve.


Fig. 14 Port connections in the LOAD and INJECT position

Sample injection

Only use a square (90°) Luer-Lock injection syringe 1/16" diameter without ridges.

The size of the syringe and the injection volume should correspond.

At least a 5-fold dose of sample solution should be used for the sample injection if working with full loop injection. This results in a relative deviation of $\leq 0.2\%$.

The sample loop can also be partially filled. Then the sample volume is determined by the injected volume.

For further details on the sample injection process see **SOP 4 Sample Injection** on page 12.

Manual 6-Port-3-Channel Switching Valves

A1367, A1368

The valves are pressure stable up to 300(400) bar. Pay attention to the general introductory sections!

Delivery Content

Order No.	Quantity	Article	Details
A1367	1	6-Port-2-Channel Injection Valve	Stainless steel, 1/16"
	6	Ferrule	1/16"
	6	Bushing UNF 10/32	Stainless steel, 12mm
	2	Open-end wrench	1/4" and 5/16"
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	
A1368	1	6-Port-2-Channel Injection Valve	PEEK, 1/16"
	3	DYNASEAL bushing	1/16" short
	3	DYNASEAL bushing	1/16" long
	6	Split-grooved clamping ring	1/16"
	6	Sealing ring	For 1/16" PEEK
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	

Layout and Function

These valves completely correspond to the other 6-Port-3-Channel Valves. The manually operated valves differentiate themselves in which they have an elongated rotor axle for feedthrough through the housing of the thermostat. In comparison to the 6-Port-2-Channel Injection Valve there is not an injection port on the 6-Port-3-Channel Valve (see Fig. 12 on page 16).

7-Port-1-Channel Valves

A1361, A1362, A1363, A1364 (manually driven)

A1373, A1374, A1375, A1376 (electrically driven)

The valves are pressure stable up to 300(400) bar. Pay attention to the general introductory sections!

Delivery Content

Order No.	Quantity	Article	Details
A1361 A1373	1	7-Port-1-Channel Switching Valve	Stainless steel, 1/16"
	7	Ferrule	1/16"
	7	Bushing UNF 10/32	Stainless steel, 12mm
	2	Open-end wrench	1/4" and 5/16"
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	
A1362 A1374	1	7-Port-1-Channel Switching Valve	PEEK, 1/16"
	3	DYNASEAL bushing	1/16" short
	4	DYNASEAL bushing	1/16" long
	7	Split-grooved clamping ring	1/16"
	7	Sealing ring	For 1/16" PEEK
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	
A1362 A1375	1	7-Port-1-Channel Switching Valve	Stainless steel, 1/8"
	7	Ferrule	1/8"
	7	Bushing UNF 10/32	Stainless steel, 12mm
	2	Open-end wrench	1/4" and 5/8"
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	
A1363 A1376	1	7-Port-1-Channel Switching Valve	PEEK, 1/8"
	7	DYNASEAL Bushing	1/8" long
	7	Split-grooved clamping ring	1/8"
	7	Sealing ring	For 1/8" PEEK
	1	DIN911 6-hexagonal key	No.3, nickel plated
	1	Manual	

Layout and Function

The four electrically operated 6-Port-3-Channel Valves are only different in material (stainless steel or PEEK) and/or in the port and a few internal dimensions. As with the 6-Port-3-Channel Valves, according to the operation mode, they are mounted to the Universal Mounting Bracket or to the Smartline Valve Drive 6 or the electrical Valve Drive K-6 (see Fig. 7 on page 11 and SOP 5 on page 14).

Fig. 15 Manual and electrical operation of the 7-Port-1-Channel Valve

The manual and electrical valves correspond in construction and functionality with each other. The essential difference is the operating lever that can be turned, clicking into place, 360°. Electrical valves, on the other hand, are moved by a Smartline Valve Drive 6 (Order No. A55020) or Valve Drive K-6 (Order No. A0377). The valves have, therefore, been described together.

All 7-Port Valves can only be switched with an Smartline Valve Drive 6 or Electrical Valve Drive K-6. Operation with a 12 or 16 positions drive results in mechanical malfunctions and may damage the HPLC system.

All 7-Port-1-Channel Switching Valves are identified by the central port 7. This port is connected with one of the six outside ports by turning it in 60° intervals.


Fig. 16 The 7-Port-1-Channel Valve switching principle

Examples for using the 7-Port-1-Channel Switching Valves

Step-function gradient

As illustrated below, the 7-Port-1-Channel Valve can be used to realise a stepwise gradient for up to 6 eluents.


Fig. 17 Stepwise gradient with 6 eluents

Fraction collector

Similarly the 7-Port-1-Channel Valve can also be used as a fraction collector in preparative chromatography (see Fig. 18).


Fig. 18 7-Port-1-Channel Valve as a fraction collector

Column switching

By combining two 7-Port-1-Channel Valves a column switching can be realised, in which up to 5 columns and a rinsing line can be involved.


Fig. 19 Column switching with two 7-Port-1-Channel Valves

In order to avoid blockage in the pump by pressure surges, it is important to always synchronously and swiftly switch the valves!

17(13)-Port-1-Channel Switching Valves

A1379, A1379-1, A1378

The valves are pressure stable up to 50(100) bar (17-port) and 25(50) bar (13-port) respectively. Pay attention to the general introductory sections!

Delivery content

Order No.	Quantity	Article	Details
A1379 A1379-1	1	17-Port-1-Channel Switching Valve	Stainless steel, 1/16"
	17	DYNASEAL Bushing	1/16" long with hexagon
	20	Split-grooved clamping ring	1/16"
	20	Sealing ring	for 1/16", injection molding polymer
	2	Open-end wrench	1/4" and 5/16"
	2	Biconical sealing ring	1/16" PETP
	1	Manual	
A1378	1	13-Port-1-Channel Switching Valve	Stainless steel, 1/8"
	13	DYNASEAL Bushing	1/8" long with hexagon
	15	Split-grooved clamping ring	1/8"
	15	Sealing ring	for 1/8" injection molding polymer
	2	Open-end wrench	1/4" and 5/16"
	2	Biconical sealing ring	1/8" PETP
	1	Manual	

Layout and Function

The 17-Port or 13-Port-1-Channel Switching Valve is operated with the Smartline Valve Drive 16 or 12 alternatively with the Electrical Valve Drive K-16 or K-12. There are no manual versions. The valves correspond, in construction and function, to the 7-Port-1-Channel Switching Valve. The outer ports are arranged in two concentric circles. Switching from one port to the next occurs in 22.5 or 30° intervals at a time.

The 17-Port Valves can only be switched with an Smartline Valve Drive 16 or Electrical Valve Drive K-16. Operation with a 6 or 12 positions drive results in mechanical malfunctions and may damage the HPLC system.

The 13-Port Valves can only be switched with an Smartline Valve Drive 12 or Electrical Valve Drive K-12. Operation with a 6 or 16 positions drive results in mechanical malfunctions and may damage the HPLC system.

Fig. 20 17-Port-1-Channel Switching Valve with Valve Drive K-16

The main area of application for the Multiport Valve is for use in fraction collection. It offers commercially available fraction collection not only at an affordable rate, but rather with functional advantages. It is essential that the size and form of the collection container is completely free to choose. Furthermore, if working with the preparative version of EuroChrom[®], the number of separated collected fractions can be extended as desired in a combination of up to 8 valves. This leaves up to 120 collection positions available.

For preparative use this means that with larger flow rates the 13-Port Valve with 1/8" connections is preferred. In this case up to 88 collection positions are available.

Fundamentally, with the 17-Port or 13-Port-1-Channel Switching Valves other functions of the 7-Port-1-Channel Switching Valve can also be fulfilled. For the production of step-function gradients or for other eluent variations this is valid without restrictions.

However, in the construction of column switching, for reliability, it is restricted to a pressure load of only 25 bar (A1378), 50 bar (A1379), or 100bar (A1379-1). The higher pressure stability of the valve, A1379-1, is achieved by a changed order of the disk springs (see Fig. 21). Changing the disk spring order the 13-Port-Valve (A1378) also can be adapted to a higher pressure stability of 50 bar. A disadvantage is that a higher amount of wear and tear must be taken into consideration.


Fig. 21 Order of the disk springs in 17-Port-1-Channel Valves

Examples for switching with Multiport Valves

In this section various examples for switching Multiport Valves are introduced. One is, however, not limited to these examples. Qualified combinations most certainly allow several other special tasks for HPLC to be realised.

Column selection with a 2-Channel Valve

The pre-separation of the sample is done on the pre-column. A part of the pre-separated sample components can be further separated on column 1 and after switching (dotted or broken lines) on column 2.


Fig. 22 Column selection with a 2-Channel Valve

Using a precolumn to precondition the sample

The sample substances are held on a precolumn. After the column switching (dotted or broken line) the sample substances are transported with another eluent or gradient in the analytical separation column.


Fig. 23 Preconditioning the sample using a precolumn and a 2-Channel Valve

Selecting a column with a 3-Channel Valve

Without interrupting the operation, the flow can be switched from column 1 onto column 2 (dotted or broken line).


Fig. 24 Selecting a column with a 3-Channel Valve

Precolumn back-flushing

Firstly the sample is guided through the precolumn before it reaches the analytical separation column. By switching (dotted or broken line) the columns are rinsed in

the opposite order and direction. This ensures that the sample components that were retained on the precolumn are eluted.


Fig. 25 Precolumn back-flushing with a 3-Channel Valve

Column back-flushing

After ending a recorded chromatogram, the column will be rinsed in the opposite direction by switching (dotted or broken line). Strongly retarded substances are removed out of the column by back-flushing with another eluent.


Fig. 26 Column back-flushing with a 3-Channel Valve

Sample enrichment

The analysis substances are enriched in the precolumn (dotted or broken line). After switching the valve, the sample is rinsed in the opposite direction of the precolumn with another eluent using a second pump. It is then transferred into the analytical separation column.


Fig. 27 Column back-flushing with a 3-Channel Valve

Alternating enrichment with a precolumn through rinsing

Pump 1 aids the sample over the injection valve to enrichment in precolumn 1. At the same time, the sample previously enriched in precolumn 2 (before the last switch) is transported with the eluent from pump 2 for analysis in the separation column. This process alternates with every switch.


Fig. 28 Alternating enrichment and precolumn cleaning with two 3-Channel Valves

In order to avoid blockage in the pump or misalignment, it is important to always synchronously and swiftly switch the valves!

Valves in which two neighboring ports are short circuited by a capillary bridge, can also be replaced with 2-Channel Valves.

Alternating enrichment with precolumn back-flushing

Strongly retarded sample components are enriched on the precolumn. After switching, these components are guided by precolumn back-flushing with another eluent to further separation on the separation column. Due to the fact of two connected precolumns it is possible without interrupting the run to proceed the enrichment alternating and parallel to the separation. For better identification the eluent current of pump 2 is represented with a broken or (after switching) dotted line.


Fig. 29 Alternating enrichment and precolumn back-flushing with two 3-Channel Valves

In order to avoid blockage in the pump or misalignment, it is important to always synchronously and swiftly switch the valves!

Valves in which two neighboring ports are short circuited by a capillary bridge can also be replaced with 2-Channel Valves.

Alternating sample preconditioning with precolumn through rinsing

The sample is preconditioned on a precolumn. Only the slightly retarded sample components are transported to the analytical column and separated there on. After switching, the precolumn is washed with another eluent in order to remove the stronger retarded components. This process alternates between two precolumns.


Fig. 30 Alternating sample preconditioning and precolumn washing with 3-Channel Valves

In order to avoid blockage in the pump or misalignment by pressure surges, it is important to always synchronously and swiftly switch the valves!

Valves in which two neighboring ports are short circuited by a capillary bridge can also be replaced with 2-Channel Valves.

Alternating sample preconditioning with precolumn back-flushing

The sample is preconditioned on a precolumn. Only the slightly retarded sample components are transported to the analytical column and separated there on. After switching, the precolumn is washed with another eluent in opposite direction in order to remove the stronger retarded components. This process alternates between two precolumns. For better identification the eluent current of pump 2 is represented with a broken or (after switching) dotted line.


Fig. 31 Alternating sample preconditioning and precolumn back-flushing with two 3-Channel Valves

In order to avoid blockage in the pump or misalignment by pressure surges, it is important to always synchronously and swiftly switch the valves.





Valves in which two neighboring ports are short circuited by a capillary bridge can also be replaced with 2-Channel Valves.

Repair and Service Note

SOP 8 Dismantling the valve body

1. Alternately loosen the three imbus mounting screws on the valve.
2. Disassemble the valve in the order shown in the exploded view illustrations (Fig. 4 and Fig. 13 on pages 10 and 17).
3. Replace the defected parts (i.e. the rotor seal).
4. Reassemble the valve in the opposite order. Pay particular attention to the correct order and alignment (i.e. disk spring).
5. The proper alignment of the rotor seal on the various valves has been illustrated in the following table.

Valve Type	Assembling position of the rotor seals
1-Channel	
Injection Valve	
2-Channel	
3-Channel	

When alternately screwing the valve together the rotor seal will automatically center itself. It is then set by the disk spring level on the stator sealing surface. Adjustment is not necessary.

Due to the special design of the valve, it is not possible to mount the rotor seal obliquely. The result is a long lifespan of the rotor seal and of the entire valve.

Rotor seals with only slight scratches can be smoothed out. Lay a piece of very fine sand paper (i.e. 1000 or rouge) on a smooth base and rub the rotor seal in even circles on it.

SOP 9 Inserting a new needle seal in the injection port

1. Remove the needle seal.
2. Push the new needle seal in the bushing.
3. Assemble the needle seal to the split-grooved clamping ring (It is customary that the old one can be reused) and a new polymer seal. Alternatively, both parts can be replaced with a polymer biconical ring.
4. Slightly tighten the bushing.
5. Insert an injection syringe into the needle seal.
6. Together with the syringe, push the the needle guide to the end of the injection port.
7. Tighten the bushing enough so that the syringe can still be easily taken out, however, with considerable friction.


Fig. 32 Injection port with needle seal

Accessories and Spare Parts

Spare Parts

Item	Valve	Order No-
Rotor seal	6-Port-valves 1/16"	A0611
	6-Port-valves 1/8"	G2527
	7-Port-valves	A0880
	17(13)-Port-valves	P3058
Reed contact	All valves	G0365
Magnetic core	All valves	M0528
Needle seal	1/16" Injection-valves	P0653
	1/8" Injection-valves	P0955

Accessories

Item	Order No-
Manual	V7447
Universal Mounting Bracket	A4018
Mounting Bracket for K-120	A0406
Analytical sample loop 1/16", stainless steel 5, 10, 20, 30, 50, 75, 100, 200, 500, 1000, 2000µl 1/16" PEEK 10, 20, 50, 200, 500µl	
Preparative sample loop 1/8", stainless steel or PEEK 1, 2, up to 11, up to 45ml 1/16", stainless steel or PEEK up to 11, up to 45ml	
Injection syringes 10, 25, 50, 100, 250, 500, 1000, 2500µl	A0723-A0730

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Guarantee statement

The guarantee period for the Switching and Injection Valves, Order No. A1357, A1358, A1359, A1360, A1361, A1362, A1363, A1364, A1365, A1366, A1367, A1368, A1369, A1369-1, A1370, A1371, A1371-1, A1372, A1373, A1374, A1375, A1376, A1377, A1378, A1379, A1379-1, and A1378 is 12 months beginning with the dispatch date from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from the guarantee. Damage caused by blockages and/or wear and tear of parts such as fuses and seals are not covered by the guarantee. Claims under this guarantee are only valid if the enclosed guarantee card has been returned to us at the address shown below within 14 days of receipt of the instrument. Defective valves should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of conformity

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**Manufacturer's name and address:
Herstellername und -adresse**

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comply with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)
EN 6100-3-2 (10/2006)
EN 61326-1 (10/2006)

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Berlin October 11, 2007

Alexander Bünz (Managing Director)

Using this Manual	5
SOPs in this manual.....	5
General	7
Capillary Connections	7
Manually driven 6-Port-3-Channel Valves	9
Delivery Content.....	9
Layout and Function.....	10
Installation	11
Operation.....	12
Using the valve as a Switching Valve	12
Electrically driven 6-Port-Multi-Channel Valves.....	13
Delivery Content.....	13
Layout and Function.....	14
Use as an Injection Valve	15
Using the valve as a Switching Valve	15
Special version for the 6-Port-2-Channel Valve.....	15
Manual 6-Port-2-Channel Injection Valve	16
Delivery Content.....	16
Assembly and Function.....	16
Manual 6-Port-3-Channel Switching Valves	18
Delivery Content.....	18
Layout and Function.....	19
7-Port-1-Channel Valves.....	19
Delivery Content.....	19
Layout and Function.....	20
Examples for using the 7-Port-1-Channel Switching Valves.....	20
Step-function gradient.....	20
Fraction collector.....	21
Column switching.....	21
17(13)-Port-1-Channel Switching Valves.....	22
Delivery content.....	22
Layout and Function.....	22
Examples for switching with Multiport Valves	23
Column selection with a 2-Channel Valve	23
Using a precolumn to precondition the sample	24
Selecting a column with a 3-Channel Valve.....	24
Precolumn back-flushing.....	24
Column back-flushing.....	25
Sample enrichment	25
Alternating enrichment with a precolumn through rinsing	26
Alternating enrichment with precolumn back-flushing.....	26
Alternating sample preconditioning with precolumn through rinsing	27
Alternating sample preconditioning with precolumn back-flushing	27
Repair and Service Note.....	28
Accessories and Spare Parts.....	30
Spare Parts	30
Accessories.....	30
Zur Benutzung des Handbuches	31
Standardarbeitsanweisungen in diesem Handbuch.....	31
Allgemeines.....	33
Kapillar-Verbindungen	34
Manuell betriebene 6-Port-3-Kanal-Ventile.....	35
Lieferumfang	35
Aufbau und Funktion	36
Installation	38
Inbetriebnahme und Bedienung.....	38
Verwendung als Injektionsventil.....	38
Verwendung als Schaltventil.....	39
Elektrisch betriebene 6-Port-3-Kanal-Ventile.....	40
Lieferumfang	40
Aufbau und Funktion	41

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Dr. Ing. Herbert KNAUER GmbH