


AZURA

▶ AZURA Compact Bio LC 50 Flow scheme

Example Configuration

- 1: Wash Buffer
- 2: Elution Buffer
- 3: Sample
- 4: CIP-Solution
- 5: Storage
- 6: Water



- 1: Waste
- 2: Flow Through
- 3: Wash
- 4: Eluted Sample
- 5: Optional
- 6: Optional

® AZURA
All rights reserved.

The information in this document is subject to change without prior notice.
Printed on environmentally friendly paper from sustainable forests.

Date of the last update: 06.07.2016
HPLC · SMB · Osmometrie

© KNAUER
Wissenschaftliche Geräte GmbH
Hegauer Weg 38
14163 Berlin, Germany

Phone: +49 30 809727-0
Telefax: +49 30 8015010
E-Mail: info@knauer.net
Internet: www.knauer.net